Guaina monouso sterile per nasofaringoscopi flessibili e rigidi (Ref. GNFxxxx; GFRxxxx; GNLxxxx)

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Attenzione

Sterile disposable sheaths for flexible and rigid nasopharyngoscopes (Ref. GNFxxxx; GFRxxxx; GNLxxxx)

Read these instructions for use carefully before each use and keep them easily accessible to interested users.

Please refer to the specific information on the intended use, indication and contraindications.

Read these instructions for use carefully. Improper use of the products may cause serious injury to the patient, user or third parties. Cover Srl rejects any responsibility for things or people in the event that other maneuvers are performed than those reported in this instruction sheet and in cases where the device was not used by qualified personnel.

USER MANUAL**DEVICE DESCRIPTION**

The COVER-ENT device is a sterile sheath to be used on the optical instrument called the nasopharyngoscope. The sheath is equipped, in the distal part, with a lens to allow the vision with the optical instrument during the clinical activity by the otolaryngologist doctor.

INTENDED USE AND CLINICAL BENEFITS

The sterile sheath for diagnostic optical instrumentation is used for the protection of optical instruments during hospital or outpatient endoscopic examinations. In particular, COVER-ENT is used in the otolaryngology field, on the instrument called nasopharyngoscope. The purpose of COVER-ENT is to act as barrier for microorganisms to prevent cross infections between patients and users due to the use of the instrument on multiple people during diagnostic examinations.

FEATURES

The cover sheath for nasopharyngoscopes COVER-ENT is available in various sizes suitable for use with the most popular nasopharyngoscope models on the market. The nominal diameters of the COVER-ENT sheath refer to intact nasopharyngoscopes and/or in any case corresponding to the specifications provided by the respective manufacturers.

The cover sheath for nasopharyngoscopes COVER-ENT is a valid barrier system to ensure simultaneously:

- The complete sterility of the nasopharyngoscope at every use
- The drastic reduction in recovery time between exams
- The reduction of procedure times

The device COVER-ENT is practical and easy to use: it fits easily on the nasopharyngoscope before use and, thanks to the reduced thickness and great elasticity of the material used, it is easily tolerated by the patient even during the orientation phases. The high transparency of the lens and its conformation ensure a clear and distortion-free vision.

The device is **STERILE** and is **DISPOSABLE**.

INDICATIONS

The use of the COVER-ENT sheath is indicated for the execution of endoscopic examinations in inpatient and/or outpatient area where nasopharyngoscopes are used, and it is reserved for adequately trained and qualified medical personnel to perform endoscopic examinations.

Can be used on both adult and pediatric patients based on the size or article code chosen.

CONTRAINDICATIONS

- Do not use the COVER-ENT sheath with damaged nasopharyngoscopes.
- The use of an unsuitable COVER-ENT sheath may damage the nasopharyngoscope.
- Do not introduce lubricants or detergents into/on the COVER-ENT sheath
- Do not force the introduction of the nasopharyngoscope into the COVER-ENT sheath. This activity may damage the sheath itself and the endoscope

PRECAUTIONS

- In the case of modified and/or repaired equipment, please check the actual diameter of the appliance before using the COVER-ENT sheath.
- Before use, check the integrity of the primary packaging (envelope). Quality and sterility can only be guaranteed for sealed and undamaged packaging. If the primary packaging is open or damaged, discard the product and use a new one.
- Check the expiry date on the label before use and do not use after the expiry date.
- Disposable device. Reuse can cause cross contamination for patient and operator.
- Do not re-sterilize. Re-sterilization may cause cross contamination for the patient and operator.
- Before use make sure you have a suitable sheath for the instrument to be covered.
- After covering the instrument with the sheath, visually verify that there are no visible cuts or holes.

INSTRUCTIONS FOR USE**PREPARATION**

Open the bag on the side of the connector, with both hands.

Check the nasopharyngoscope and make sure it is clean and dry, otherwise, gently wipe the surface with a detergent and dry thoroughly.

Gently introduce the distal end of the nasopharyngoscope into the sheath.

POSITIONING

Slide the sheath over the instrument and, through a slight pull, attach the connector to the conical part of the nasopharyngoscope.

Remove the protective tube and store it until the end of the examination

REMOVAL

Reposition the protective tube and remove the sheath by slightly rotating the connector. Remove the sheath from the distal portion of the nasopharyngoscope. In case of resistance, the sheath can be removed by pulling it from the part of the lens. If the sheath is glued, place a small amount of alcohol with a needle-free syringe inside the sheath. After removal, dispose of the sheath with the protective tube.

REPORTING OF MAJOR ACCIDENTS

The user and/or the patient is recommended to report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

CONSERVATION

If stored in the original packaging, the medical device is valid for 5 years. The medical device can be stored and transported in a dry place, in its original packaging. Also protect from sunlight, heat and handle with care.

DISPOSAL METHODS

Disposable medical device. After use, throw the product into the appropriate hospital containers and according to the protocols of the health facility where the product is used. Do not disperse in the environment.

SYMBOLS USED ON THE LABELS OF THE PRODUCT

Do not reuse- Disposable

Medical Device

Expiry date

Do not use if the envelop is open or damaged

Keep away from light and heat sources

Manufacturer

Date of manufacture

Sterilized by EtO. Single layer sterile barrier system

Refer to the instructions for use

Keep dry

Product code

Lot number

Do not re-sterilize

Medical device according to the Regulations UE 2017/745

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